

Product Name :
Meat Training Equipment

Product Code :
NLAB-EDUCATIONTVET29001



Description :

Meat Training Equipment

Technical Specification :

Device adapt the stainless steel and Al-Mg material, strong and durable.

Small body, easy to cleaning, convenient to carry.

Dozens of blade combination meat shortcut.

Input power supply: two phase three wire system 220V±10% 50Hz

Boundary dimension: 600mm×400mm×500mm

Motor power? 0.5KW

Weight??30kg

Working condition: environment temperature: -10??+40

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003