

**Product Name :**  
Mechanical Stop Clock

**Product Code :**  
NLAB-EDUSCHOOLEQ46002



**Description :**

Mechanical Stop Clock

**Technical Specification :**

A stop clock is a small portable device which used to measure the amount of time between the start and end of a particular task.

It has a start and stop lever.

The dial has a second's and minute's pointer to indicate the time period in the most accurate manner.

The clock is also fitted with a fly back lever better reading and accuracy

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003