

| | |
|---|--|
| Product Name : Medical Compressor | Product Code : NLABMEDICAL-M32100011 |
|---|--|

| |
|--|
| |
|--|

Description :

Medical Compressor India, Medical Compressor manufacturer India, suppliers, Medical Medical Compressor India, Hospital Medical Compressor India

Technical Specification :

Medical Compressor
H.P: 0.5/0.9
Volt: 230V
Tank: 38 litre
Weight: 36kg
ABAC (oil free)
H.P: 2
Volt: 230V
Tank: 50 litre
Weight: 38kg

Naugralabequipments
Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003