

**Product Name :**  
Medical Infrared Thermometer

**Product Code :**  
NLAB-SAFETYLABP11016



**Description :**

Medical Infrared Thermometer

**Technical Specification :**

With Celsius/Fahrenheit switch module

Non-contact measurement, safe without exposure to radiation

Ergonomic design for comfort grip

<b>Measurement Range:</b>
Normal (Body Core) $\pm 0.3^{\circ}\text{C}$ (32-34.9 $^{\circ}\text{C}$ and 42.1-42.2 $^{\circ}\text{C}$ )
<b>Resolution:</b>
<b>Memory Function Storage:</b>
<b>Power Supply:</b> (DC3V)
<b>Dimensions:</b> 87x138mm
<b>Weight:</b> 37g (without Batteries)
<b>Automatic Shutdown:</b> 1min

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003