

Product Name :
Medical Syringe Pumps

Product Code :
NLABMEDICAL-M3420002

Description :

Syringe Pumps India, Syringe Pump manufacturer India, suppliers, Medical Syringe Pump India, Hospital Syringe Pump India

Technical Specification :

Freely stackable: Anyone freely stack one onto another to provide multiple solutions.

Combination solutions that have a wide range of clinical applications.

Three Power Source: AC22050Hz, Internal Battery, 12V car charge.

Digital LCD display: High capacity words, friendly user interface, dynamically display working status.

Positive and Negative Shuttle Debugger:

Can freely regulate parameters.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003