

Product Name : Medical Syringe Pumps	Product Code : NLABMEDICAL-M3420002
Description : Syringe Pumps India, Syringe Pump manufacturer India, suppliers, Medical Syringe Pump India, Hospital Syringe Pump India	
Technical Specification : Freely stackable: Anyone freely stack one onto another to provide multiple solutions. Combination solutions that have a wide range of clinical applications. Three Power Source: AC22050Hz, Internal Battery, 12V car charge. Digital LCD display: High capacity words, friendly user interface, dynamically display working status. Positive and Negative Shuttle Debugger: Can freely regulate parameters.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003