


Product Name :
Megger

Product Code :
NLAB-TESTINGLABM19003



Description :

Megger

Technical Specification :

Technical Specifications :-

Range : 200M Ω , 100V, 200M Ω , 250V, 200M Ω , 500V, 1000M Ω , 1000V

Accuracy : $\pm(3\% + 1d)$; Resolution : 0.1M Ω , 1M Ω

Sampling time : 0.4 sec.

Display : Minimum 0.7" LCD & Max. Indication : 1999

Measurement with high drive current, short circuit current > 2.55mA

LSI-circuit for high reliability and durability

Built in over input and low battery indication

Overload protection for each range

Automatic zero adjust and automatic circuit discharge

Portable housing case with front protective cover

Power : DC 9 V battetry.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003