

**Product Name :**  
Membrane Filter Holder – All Glass

**Product Code :**  
NLAB-SCIENTIFICEQ18007



**Description :**

Membrane Filter Holder – All Glass

**Technical Specification :**

**Borosilicate Glass**, with Sintered disc for membrane support and aluminium 'duck' clamp.

- Membrane Filter Holder – All Glass. For Membrane diameter 47mm. Funnel Capacity 300 ml, Flask Capacity 1000 ml.
- Membrane Filter Holder – All Glass. For Membrane diameter 25mm. Funnel Capacity 35 ml, Flask Capacity 250 ml.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003