

Product Name :
Metacentric Height Apparatus

Product Code :
DM198



Description :

Metacentric Height Apparatus

Technical Specification :

Fundamentals of Buoyancy, Metacenter and Metacentric Peak are widely utilized in Design and steadiness of Ships, Submarines, Scorching air balloons and so forth. Equipment encompass SS reservoir tank set up on sturdy body. A small ship adaptation made from Perspex has a provision to position load on the deck and has an association to measure the inclination as a result of loading. Optionally test may also be additional modified to characterize stipulations as cargo ship and battle ship. Specifications: Reservoir Tank Measurement: 500 X 500 X 200 mm Ship Adaptation Dimension: 230 X 350 mm Round Weights: four No's Graduated Arc for measuring tilt perspective Specific Technical Guide and On-web page Coaching

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003