

Product Name :
Metallurgical Samples

Product Code :
NLAB-MECHANICALAB250016



Description :

Metallurgical Samples

Technical Specification :

Bakelite Mounted Ready Specimen for microscopic Observation, A set consists of 23 Standard samples which covers the entire range. Complete information including chemical composition, mechanical treatment, details of heat treatment, type of Etchant with description of Microstructure along with their colored photographs are provided in a booklet.

FEATURES

These sets are supplied in a wooden storing box with silica gel.

Type of Etchant with description of Microstructure along with their colored photographs are provided in a booklet.

LIST OF METALLIC MICROSTRUCTURE

Dead Mild Steel

Mottled cast Iron

Low Carbon steel
Ductile cast Iron
Medium Carbon steel (Annealing)
Cartridge Brass
Medium Carbon steel (Normalizing)
Muntz Metal
Decarburized High Carbon steel
Tin Bronze
Inclusion in steel
Electroplated Component
Hardened steel
Anodized Aluminum
Tempered steel
Fusion Welded Mild Steel
Carburized steel
Friction Welded Steel
Tool steel
Powder Metallurgy Component
Gray cast Iron
Deformed Mild Steel
White cast Iron

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003