

Product Name :
Alcohol Thermometer

Product Code :
NLAB-SCIENTIFICEQ14001



Description :

Alcohol Thermometer

Technical Specification :

Mercury in glass capillary, yellow backed with scale printed in black for easy reading. 5-7mm diameter with reinforced bulbs. Packed in cardboard or plastic case. Range Graduation Approx Length

- 0 to 50°C 1° 300mm
- -10 to 110°C 1° 300mm
- 0 to 150°C 1° 300mm
- 0 to 200°C 1° 300mm
- 0 to 250°C 1° 300mm
- 0 to 300°C 1° 300mm
- 0 to 360°C 2° 300mm

- 0 to 400°C 2° 300mm
- 0 to 220°C 1° 300mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003