

<b>Product Name :</b> Microglass Fiber Filter Paper	<b>Product Code :</b> NLAB-CHEMICALSIND19008
<b>Description :</b>	
Microglass Fiber Filter Paper	
<b>Technical Specification :</b>	
<p>Glass Micro Fiber Filters manufactured from borosilicate glass. Because of their intrinsic properties, Glass Micro Fiber Filters have wide applications in many areas of Laboratory Analysis, especially when fine filtration and high loading capacity is required. Microglass fiber filters are also used in the filtration of hot gases and liquids, and in pre-filtration. Applications for Glass Micro Fiber Filters include environmental analysis for water, wastewater, testing of soils, air pollution and monitoring, research and process biochemistry, and gravimetric analysis involving ignition of samples.</p>	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003