

Product Name :
Microprocessor Based Automatic Tissue Processor

Product Code :
NLAB-LABORATORYEM28019



Description :

Microprocessor Based Automatic Tissue Processor

Technical Specification :

Feather touch keypad.

Digital display for adjusting time sequence.

Battery backup to retain data.

12 stage timing sequence can be programmed with each step of 1 minute to 16 hours 39 minutes (999 min.)

Delay action of 1 hour to 99 hours.

2 Electronically controlled wax bath 10 deg C to 80 deg C

1 or 2 Stainless Steel Basket.

24 tissue containers (Cassettes).

1 or 2 stainless steel basket rotors.

10 without lip glass beakers ($\frac{1}{2}$ ltr, 1 ltr or 2 ltr).

Accessories (Optional) :-

Additional basket rotor.

Tissue capsule basket SS.

Tissue cassettes with sliding cover.

S. Shape divider for tissue capsule.

V. Shape divider for tissue capsule.

Additional wax bath.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003