

**Product Name :**  
Microprocessor Controlled Slide Staining Machine

**Product Code :**  
NLAB-LABORATORYEM28020



**Description :**

Microprocessor Controlled Slide Staining Machine

**Technical Specification :**

Feather touch keypad.

Digital display for adjusting timing sequence.

Battery backup to retain data.

12 stage timing sequence can be programmed with each step of 1 minute to 16 hours 39 minutes (999 min.)

12 Staining troughs

Compact bench top model.

Perfect H & E stainer for routine application.

High specimen through put.

Minimum reagent exposure.

Easy to clean.

Simple to operate.

Highly recommended for hematology and also for staining of blood films and bone marrow spreads etc.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003