

**Product Name :**  
Microprocessor Controlled Tissue Processor

**Product Code :**  
NLAB-LABORATORYEM28021



**Description :**

Microprocessor Controlled Tissue Processor

**Technical Specification :**

Feather touch key pad

Digital display for adjusting time sequence

Battery backup to retain data

Stage timing sequence can be programmed with each step of 1 minute to 16 hours 39 minutes (999 min.)

Delay action of 1 hour to 99 hours

2 Electronically Controlled wax bath 10°C to 80°C

1 or 2 Stainless Steel Basket

24 tissue containers (Cassettes)

1 or 2 Stainless steel Basket Rotors

10 without lip glass beakers (½ ltr, 1 ltr or 2 ltr.)

Accessories :

Additional Basket Rotor

Tissue Capsule Basket SS

Tissue Capsule with Sliding Cover

S Shape divider for Tissue Capsule

V Shape divider for Tissue Capsule

Additional Wax bath

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003