

**Product Name :**  
Microprocessor pH System

**Product Code :**  
NLAB-LABORATORYEM15001



**Description :**

Microprocessor pH System

**Technical Specification :**

Measure pH, mV & Temperature

Highly Accurate & Easy to Operate

Auto Temperature Compensation Facility

7 Segment Bright Red LED Display

Soft touch membrane keys

Microprocessor pH System is a solid state instrument designed to provide the precise pH, mV and temperature measurements.

It has auto buffer recognition for 4.00, 7.00 and 9.2 pH buffers and auto 2 or 3 points buffer calibration.

It has 7 segment bright red LED display.

The instrument uses the latest microprocessor technology and advanced engineering techniques so as to give enhanced accuracy and reproducibility.

The system has 6 soft touch membrane type keys for ease of operation.

The calibration is retained even after the instrument is switched OFF.

The instrument is extremely useful for agricultural and soil analysis laboratories, swimming pools, water quality control in boiler feed water, water works department, fertilizer plants, petroleum refineries, breweries, water purification plants etc.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003