


Product Name : Microscope Eyepiece Camera	Product Code : NLAB-EDUSCHOOLEQ173014
---	---



Description :

Microscope Eyepiece Camera

Technical Specification :

Works with both monocular and binocular microscopes by replacing the eyepiece.

Basic microscope eyepiece camera with 2 MP resolution, which connects to a PC via the inbuilt USB 2.0 cable.

Fits eyepiece tubes with a diameter of 23 mm or 30 mm (adapter included).

Comes with software for Windows PCs.

640 x 480 live resolution

2 MP still resolution

Naugralabequipments
Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003