

<b>Product Name :</b> Microscope Eyepiece Camera	<b>Product Code :</b> NLAB-EDUSCHOOLEQ173014
<b>Description :</b>	
Microscope Eyepiece Camera	
<b>Technical Specification :</b>	
Works with both monocular and binocular microscopes by replacing the eyepiece.	
Basic microscope eyepiece camera with 2 MP resolution, which connects to a PC via the inbuilt USB 2.0 cable.	
Fits eyepiece tubes with a diameter of 23 mm or 30 mm (adapter included).	
Comes with software for Windows PCs.	
640 x 480 live resolution	
2 MP still resolution	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003