

Product Name :
Microtomes Lab Semi Automatic Rotary Microtome

Product Code :
NLAB-EDUSCHOOLEQ29023



Description :

Microtomes Lab Semi Automatic Rotary Microtome

Technical Specification :

It comes with an advanced imported drive system for smooth trimming and sectioning, LCD screen displays section and trimming thickness range, number of sections etc.

It comes with a lockable hand wheel in any position, huge volume waste tray for easy operation, interchangeable specimen clamp along with a safety alarm system.

Its high precision calibrated controls make accurate positioning for re-cuts possible.

Technical Specifications Section Thickness Setting Range 0.25~60 & micro's Trimming Thickness Setting Range 0.25~60 & micro's Setting values 0.25~ 1- μ 's increment step 0.25 μ 's 1~10 micro's increment step 1 μ 's 10~ 20- μ 's increment step 2 μ 's 20~ 60- μ 's increment step 5 μ 's Minimum setting value 0.25 μ .

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

