


Product Name :
Microwave Test Bench

Product Code :
NLAB-ELECTRONICSAB200004



Description :

Microwave Test Bench

Technical Specification :

A selection of wave-guide components

The power supply for the microwave source

A detector

A meter, which monitors the detector output

Wall Chart

Microwave Chart

These Training Benches are completely self-contained & provide the means to allow students to carry out practical work at extremely low cost. A comprehensive manual containing extensive microwave theory and a progressive series of assignments is supplied with the Trainer.

Experiments:

Study of E-Plane, H-Plane & Magic Tees

Study of Directional Couplers

Study of Fixed and Variable Attenuators

Study of Isolator and Circulators

Instruments:

SKPS----- 1

VSWR Meter ----- 1

Components:

Cooling Fan----- 1

C.D.Coupler 20 dB -----1

Detector Mount----- -1

D.R.F.Meter----- 1

E-Plane Bend -----1

E-Plane Tee ----- -1

Fixed Attenuator 10 dB -----1

Fixed Attenuator 3 dB -----1

Fixed Attenuator 6 dB -----1

H-Plane Bend----- 1

H-Plane Tee ----- -1

Isolator----- 1

Klystron Mount ----- 1

Magic Tee ----- 1

Matched Termination ----- -2

Movable Short-----	1
M H D Coupler 10 dB-----	1
M H D Coupler 3 dB -----	-1
S S Tuner-----	1
Slotted Section -----	1
T Circulator -----	1
Y Circulator-----	1
Tunable Probe -----	1
Variable Attenuator 20dB -----	1
Wave Guide Stand -----	-3

Accessories:

Mains Cord -----	2
Cable BNC - BNC -----	-2
Op. Manual -----	-1

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003