

Product Name :
Mini Centrifuge 5

Product Code :
NLAB-EDUSCHOOLEQ45026



Description :

Mini Centrifuge 5

Technical Specification :

5 microcentrifuge includes a four position rotor for compatibility with a variety of 5ml test tubes, including 5ml snap cap centrifuge tubes.

The microtube positions are designed to accept 0.5 and 0.2ml tubes with the use of optional adapters.

Open the lid, and the rotor decelerates to a stop for removing samples.

In addition, this combination rotor has four additional slots for centrifuging standard microtubes (1.5/2.0ml).

The 5ml rotor positions accept 5ml centrifuge tubes (snap cap) directly.

Simply close the lid and the centrifuge quickly accelerates to the fixed speed of 5,500rpm.

While 5ml blood and culture tubes (12x75mm) and cryovials are also accepted with the included adapters.

Speed: 500 rpm / 2,000 xg	
Capacity: 4 x 1.5/2.0ml	
Dimensions (W x D x H): 14 x 20 x 11.2 cm	
Weight: 1 kg	
Electrical: Voltage 100-240V, 50-60 Hz, 300W	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003