

Product Name :
Mini Centrifuge

Product Code :
NLAB-EDUSCHOOLEQ45031



Description :

Mini Centrifuge

Technical Specification :

The mini centrifuge is a personal centrifuge that virtually fits in the palm of your hand.

Simply close the lid and the rotor instantly reaches 6,000 rpm.

When the rotor is not in use, storage in this compartment prevents loss or misplacement of the rotor.

Open the lid, and the rotor quickly decelerates to a stop for removing samples.

An increased capacity, eight-position microtube rotor is included and comes installed.

Ideal for quick spin downs of microtubes and PCR tubes, the is extremely easy to use.

Also included is a second rotor for 0.2ml PCR tubes and strips which is stored (together with the rotor removal key) in a convenient compartment on the bottom of the centrifuge.

6,000 rpm / 2,000 xg	
Cap: 50ml tubes	
32 x 0.2ml tubes	
4 x PCR Strips	
Dimensions (W x D x H) x 15 x 11.4 cm	
Weight: 1 kg	
Electrical Voltage	
100-240V, 50-60 Hz, 250W	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003