

Product Name :
Miniature Viewing Box

Product Code :
NLAB-EDUSCHOOLEQ45041



Description :

Miniature Viewing Box

Technical Specification :

This economical lightbox is suitable for viewing two stained gels.

This unit can be operated with 6 "AAA" batteries or with a 9V adapter.

Using the latest cold cathode technology makes it incredibly lightweight.

The unit features a durable plastic case, On/Off switch, and a stain resistant 4" x 5" viewing surface.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003