

**Product Name :**  
Mobile Examination Operation Theatre Lamp

**Product Code :**  
NLABMEDICAL-M2750004

**Description :**

Mobile Examination Operation Theatre Lamp

**Technical Specification :**

Dome Head: 504mm

Halogen Bulbs: 24V x 150 W x 2 nos.

Light Intensity at 1 mt: 90,000 Lux

Height Adjustment: 1078 mm

Dichoric Glass Reflector: 390mm

Color Temperature: 4200 K  $\pm$ 300

Temp rise in field: 2°-3°C

Control Unit SMPS: Available with Battery Backup

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003