

**Product Name :**  
Mobile Led Operation Theatre Lamp

**Product Code :**  
NLABMEDICAL-M2750003

**Description :**

Mobile Led Operation Theatre Lamp

**Technical Specification :**

Easy to move on low inertia, anti static, castors. The light head is mounted on a spring – balanced arm for easy position.

**Emergency Battery Back up**

Located in the base, a sophisticated battery provides up to 3 hours of uninterrupted run-time in the event of power loss. The Batteries are equipped with a “quick charge” facility allowing ultra fast battery charging.

**Extended Arm Reach**

The unique design of the Mobile base allows the light to be manoeuvred much closer to the table without impeding the surgeon. Small and compact, the reduced footprint still offers complete stability.

**Technical Specification:**

Illumination Intensity at 1m: 160000

Field size of 1m distance: 17 – 25 cm

Color Temperature: 4300 K

Working Distance: 120 cm

Color Rendering Index: 93

Number of Light Modules: 28

Dome Head Dia: 60 cm

LED Lifetime: 50000 h

Light Rotation: 360°

Control Unit: SMPS

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003