

<b>Product Name :</b> Modular Chapati Maker	<b>Product Code :</b> NLAB-VOCATIONALTVET27045
--	---

--

**Description :**

Modular Chapati Maker

**Technical Specification :**

Material: Stainless Steel
Automation Grade: Automatic
Frequency: 60 Hz
Voltage: 240 V
Power: 2.5 Kw

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

