

|  |  |
|--|--|
| <b>Product Name :</b><br>Modular Evaporator Family   | <b>Product Code :</b><br>NLAB-TECHNICALAB37008 |
|  |  |
| <b>Description :</b>   |  |
| Modular Evaporator Family  |  |
| <b>Technical Specification :</b>   |  |
| A laboratory evaporation system capable of being configured as rising or falling film, single or double effect.  |  |
| Temperature controlled pre-heat stage.   |  |
| Adjustable re-circulation on each evaporation stage.   |  |
| Manual control console.  |  |
| Integral USB interface for computer data logging and control.  |  |
| Contains vacuum pump, condenser and condensate vessel.   |  |
| Options to operate from external steam supply, or integral electrically heated pressurised hot water circulator. |  |

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003