

Product Name : Modular Evaporator Family	Product Code : NLAB-TECHNICALAB37008
--	--

--

Description :

Modular Evaporator Family

Technical Specification :

A laboratory evaporation system capable of being configured as rising or falling film, single or double effect.

Temperature controlled pre-heat stage.

Adjustable re-circulation on each evaporation stage.

Manual control console.

Integral USB interface for computer data logging and control.

Contains vacuum pump, condenser and condensate vessel.

Options to operate from external steam supply, or integral electrically heated pressurised hot water circulator.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003