

Product Name :
Monitor Trolley

Product Code :
NLABMEDICAL-M2620007



Description :

Monitor Trolley

Technical Specification :

Specification :-

Monitor Trolley is made of Mild Steel Sheet.

It has one drawer & 4 shelves.

It comes with oxygen cylinder cage & adjustable I.V. Rod.

The Monitor Trolley is mounted on 12.5cm castors.

It has a switch & socket panel for medical equipment.

Product Size : 80 L x 60 D x 140 H cms.

We are Manufacturer, Exporters and Suppliers of Monitor Trolley at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003