

Product Name :
Monochrome Display Finger Tip Oximeter

Product Code :
NLABMEDICAL-M3410002

Description :

Monochrome Display Finger Tip Oximeter India, Monochrome Display Finger Tip Oximeter manufacturer India, suppliers, Medical Monochrome Display Finger Tip Oximeter India, Hospital Monochrome Display Finger Tip Oximeter India

Technical Specification :

Color Display Finger Tip Oximeter
High quality 65k color OLED display
Automatic power-off without signal for 5 seconds
Superior accuracy and durability
Two "AAA" Alkali battery and OLED bar graph
Two display direction adjustable

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003