

**Product Name :**  
Monocular Coaxial Focusing Beginner Microscope

**Product Code :**  
NLAB-EDUSCHOOLEQ1730116

**Description :**

Monocular Coaxial Focusing Beginner Microscope

**Technical Specification :**

The Microscope offers great value and performance and features a durable cast metal frame and LED illumination.

They are achromatic, color-coded, parcentered and parfocaled, which offers a sharp clear image.

Monocular microscope with 45° inclined head that can rotate 360°.

Coaxial coarse and fine focus controls with slip clutch to prevent damage to gears

6 Hole disc diaphragm.

Locked on 10X wide field eyepiece with pointer.

These beginner microscopes are commonly used in schools for grades 4 to 8 and are excellent for young children who are studying biology and perfect for beginners. Combining with a 10X eyepiece, the magnification is 40X, 100X and 400XR.

Reverse nosepiece with 3 achromatic, parfocal and parcentered objectives: 4X, 10X, 40XR (retractable).

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003