

Product Name : Monocular Inclined Microscope	Product Code : NLAB-EDUSCHOOLEQ1730162
--	--

--

Description :

Monocular Inclined Microscope

Technical Specification :

Condenser: Abbe Condenser N.A.1.25, with iris Diaphragm with swing out filter holder.

Interpupillary Distance: 50-75mm.

Viewing Head: Monocular Inclined at 45°, Rotatable 360°.

Objective: Color coded DIN Objectives- 4x,/0.10, 10x,/0.25, 40x,/0.65 (S), 100x,/1.25 (Oil).

Illumination: Built in LED Cooled Bright Illumination or 6V 20W.

Focusing: Separate Coarse & fine Focusing Graduated knobs.

Working Stage: Large 120x140 mm Built-in low position coaxial mechanical stage.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003