

Product Name : Monocular Inclined Microscope	Product Code : NLAB-EDUSCHOOLEQ1730162
Description :	
Monocular Inclined Microscope	
Technical Specification :	
Condenser: Abbe Condenser N.A.1.25, with iris Diaphragm with swing out filter holder.	
Interpupillary Distance: 50-75mm.	
Viewing Head: Monocular Inclined at 45°, Rotatable 360°.	
Objective: Color coded DIN Objectives- 4x,/0.10, 10x,/0.25, 40x,/0.65 (S), 100x,/1.25 (Oil).	
Illumination: Built in LED Cooled Bright Illumination or 6V 20W.	
Focusing: Separate Coarse & fine Focusing Graduated knobs.	
Working Stage: Large 120x140 mm Built-in low position coaxial mechanical stage.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003