

Product Name : Monocular Laboratory Microscope	Product Code : NLAB-EDUSCHOOLEQ1730150
Description :	
Monocular Laboratory Microscope	
Technical Specification :	
<p>The precision coaxial focusing and low-position coaxial stage controls are a pleasure to use and deliver smooth, reliable performance.</p> <p>The higher quality semi-plan DIN optics provide superb optical clarity at 40, 100, 400, and 1000x magnification with 80-85% flatness of field.</p> <p>The bright 3-watt LED lighting with variable intensity control, 1.25 N.A. Abbe condenser, and iris diaphragm with swing out filter holder provide great contrast control.</p> <p>Other features like the reverse position objective turret facilitate slide changing and ease of operation.</p>	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003