


Product Name : Monocular Laboratory Microscope	Product Code : NLAB-EDUSCHOOLEQ1730150
--	--



Description :

Monocular Laboratory Microscope

Technical Specification :

The precision coaxial focusing and low-position coaxial stage controls are a pleasure to use and deliver smooth, reliable performance.

The higher quality semi-plan DIN optics provide superb optical clarity at 40, 100, 400, and 1000x magnification with 80-85% flatness of field.

The bright 3-watt LED lighting with variable intensity control, 1.25 N.A. Abbe condenser, and iris diaphragm with swing out filter holder provide great contrast control.

Other features like the reverse position objective turret facilitate slide changing and ease of operation.

Naugralabequipments
Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003