

<b>Product Name :</b> Monocular Upright Metallurgical Microscope	<b>Product Code :</b> NLAB-EDUSCHOOLEQ1730141
<b>Description :</b>	
Monocular Upright Metallurgical Microscope	
<b>Technical Specification :</b>	
Illumination: Brightfield vertical illuminator 6V-20W Halogen with field diaphragm, filter slots and bulb centring mechanism.	
Objectives: Achromatic DIN M4X-M10X & M40X.	
Stage: Large Co-axial Mechanical Stage for X-Y movement of specimens.	
Eyepieces: H6X-WF10X (any two)	
Focusing: Both side co-axial focussing knobs.	
Body: Monocular Head inclined at 45-degrees.	
Nosepiece: Quadruple revolving nosepiece with accurate centering & positive click stops.	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003