

**Product Name :**  
Mortise Machine

**Product Code :**  
NLAB-MACHINELABWK33009



**Description :**

Mortise Machine

**Technical Specification :**

Industrial

Rotates through 180o door lock morticing

Fully adjustable to sideways

Floor Standing

Power: Input 1300W

Chisel stroke: 140mm

Centre of chisel to back fence: 100mm

Max. height of timber with 12.7mm chisel and bit: 190mm

Max chisel size (softwood): 25mm

Max. chisel size (hardwood): 19mm

Type: Floor Standing

Must come complete with bushes for ¾", 13/16" and 1.13/16" shank chisels; cabinet base; service tools.

Wood Workshop Portable Tools Mortise Machine Manufacturer, Supplier and Mortise Machine Exporter in India

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003