

Product Name :
Motorcycle Electronic Fuel Injection System Trainer

Product Code :
TK188



Description :

Motorcycle Electronic Fuel Injection System Trainer

Technical Specification :

Motorcycle Electronic Fuel Injection System Trainer

The Motorcycle Electronic Fuel Injection System Trainer includes a complete OEM motorcycle electronic fuel injection system and engine management. It can be based on any model of your choice and designed to incorporate any components you require.

Educational Advantages of Motorcycle Electronic Fuel Injection System Trainer

Injection system including four (4) injectors on OEM rails

Electronic control module (ECM)

Data Link Connector (DLC) for complete diagnosis

120-V AC/12-V DC supply unit with protective circuit

DC electric motor to simulate the engine speed (RPM)

CKP and CMP assemblies

Fan motor

Speedometer

Control panel

OEM fuse block

Eight (8) liters of injection liquid

Injected liquid pan

Fuel reservoir

With Full Accessories.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003