


Product Name :
Motorised Casangarnde Equipment

Product Code :
NLAB-VOCATIONALTVET18005



Description :

Motorised Casangarnde Equipment

Technical Specification :

Gear Motor: Attached to the metal base; enclosed in the housing with the on-off switch.

Counter: Mechanical; records number of cup drops.

Cup. Spun brass; pin and slot mounting.

Base. Molded hard rubber.

Cam. Molded, self-lubricating plastic.

Housing. Aluminum; holds the cup and drop adjustment parts.

Power requirements 220v AC, 50 Hz.

Grooving Tools.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003