

Product Name :
Motorised CBR Machine

Product Code :
NLAB-VOCATIONALTVET18016



Description :

Motorised CBR Machine

Technical Specification :

CBR test for determining bearing capacity of soils.

Consisting twin column frame and a motorized drive system.

Capacity: 50 kN.

Load Speeds: Selectable: 1.27 mm (0.5") or 1.0 mm.

Horizontal Clearance: 255 mm (10").

Vertical Clearance: 800 mm (31.5") maximum.

Platen size: 133 mm (5.24") diam.

Platen Travel: 105 mm (4.13")

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003