

**Product Name :**  
Motorised Rockwell System Hardness Tester

**Product Code :**  
NLAB-TESTINGLABM13004



**Description :**

Motorised Rockwell System Hardness Tester

**Technical Specification :**

- Suitable for Rockwell Tests.
- Motorised for automatic Operation cycle i.e. load/dwell/unload.
- Preliminary Test Force - 98.07N (10kgf).
- Additional Test Force - 490.3, 882.6, 1373N (50, 90, 140 kgf).
- Total Test Force - 588.4, 980.7, 1471N (60, 100, 150 kgf).
- Test Force Selection by external dialing.
- Auto zero setting dial gauge.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003