

**Product Name :**  
Motorized Chain Pulley Block

**Product Code :**  
NLAB-MECHANICALAB240043



**Description :**

Motorized Chain Pulley Block

**Technical Specification :**

0.5 Ton to 5.0 Ton with GR-80 alloy steel load chain

Available with single phase, three phase motors

With mechanical, electromagnetic brake

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003