

Product Name :
Motorized Compression/Tension Machine

Product Code :
NLAB-VOCATIONALTVET17001



Description :

Motorized Compression/Tension Machine

Technical Specification :

Use:-

For Compression & Tension Tests Machine For Tensile Test (500Kn) And Compression Test (2000 Kn).

Max. Tensile Test Load: 500 Kn (0.01 Kn Divisions)

Max. Compression Test Load: 2000 Kn (0.1 Kn Divisions)

Piston Stroke: 100 Mm

Span Between the Columns: 280 Mm

Manual Regulation of the Load/Tensile Rate With Pacer

Accessories:-

Compression Load Cell 2000 Kn Capacity (En 2) With Factory Calibration Certificate.

Two Channel Universal Digital Indicator For Dynamometers And Load Cells.

Remote Button Switch For Readout Units.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003