

Product Name : Motorized Research Microscope	Product Code : NLAB-EDUSCHOOLEQ173057
Description : Motorized Research Microscope	
Technical Specification : Condenser: Motorized swing-out abbe, universal Five/Six position Nosepiece: Motorized Septuple/Sextuple nosepiece(Auto Koehler's Illumination) with slots for POL / DIC Eyepiece: Super wide field Optical system: Universal Infinite optical system, color corrected Viewing Head: Inclined, Two/Three position ergo head Koehler's Illumination: S-LED, Koehler Illumination, Eco-Energy Saving, Auto Cut Off, Fly Eye Lens built-in filters Mechanical Stage: Motorized Double Layer, Rackless Ceramic Hard Coated, height adjustable co-axial movemet knob, Stage Lock, Dual Slide Holder, single hand operation Objectives: Plan, Phase Achromatic, Fluotar, SEMI-APO, APO Darkfield Attachment: Cardioid oil condenser all magnification range Fluorescence Lamp: High intensity mercury lamp, cool LED, fiber optic illumination Fluorescence Filter: Wide, Narrow Band Pass, Multi Band DIC: Slot for holding DIC plate and prism Teaching Head: Upto 10 person multi viewing attachment Filters: Built-in ND Filters Motorized: Nosepiece, focuse, condenser, fluorescence Focusing: 3/2 step focusing	

Fluorescence: Six filter turret holder

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003