

**Product Name :**  
Multi-Channel Temperature Data Logger With 16 Channels

**Product Code :**  
DM993



**Description :**

Multi-Channel Temperature Data Logger With 16 Channels

**Technical Specification :**

Multi-Channel Temperature Data Logger With 16 Channels

A multi-channel temperature recorder controlled by ARM microprocessor.

Color display and multi-channel parallel test.

Multi-channel temperature recorder can collect, alarm and transmit multi-channel temperature simultaneously.

Wide adaptability and supporting K/N/E/J/T/R/S/B thermocouple.

Test range is from -200<sup>o</sup>C to 1800<sup>o</sup>C.

Multi-channel temperature recorder equipped with interface that can directly upload data to PC, and also has USB interface.

When it is not convenient to connect to PC, it is save measured data directly to USB disk, and then transfer data to PC when needed.

Accept module combination.

Data acquisition, analysis and printing is easily realized through PC software.

With a strong anti-vibration shell.

The USB interface can store the collected data in the U disk through the computer monitoring software

Include Batteries, USB Cable, PC Software, Card Slot Cover and User Manual.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003