

**Product Name :**  
Multi Functional Patient Care Nursing Manikin

**Product Code :**  
NLAB-MEDICALMODELAB11006



**Description :**

Multi Functional Patient Care Nursing Manikin

**Technical Specification :**

The model is made of soft and half-hard PVC plastic, solid and durable, easy to use.

All joints are movable, can be bent, all parts.

The model is made of soft and half-hard PVC plastic, solid and durable, easy to use.

All joints are movable, can be bent, all parts can be disassembled.

It can demonstrate following items:-

Washing face and bathing the bed-ridden patient

Care and cleaning of the mouth

Care and cleaning of trachea dissection

Inhalation of oxygen (nasal cork or nasal catheter)

Nasal lavage

Stomach lavage

Heart infection

Heart resuscitating on chest

Chest puncture

Liver puncture

Kidney puncture

Abdomen puncture

Medulla puncture

Lumbar puncture

Deltoid injection

Deltoid hypodermic injection

Vein injection

Vein puncture

Vein transfusion

Vein blood transfusion

Urethral catheterization of women

Muscular injection on buttocks.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003