

Product Name :
My Block Mini Digital Dry Bath

Product Code :
NLAB-EDUSCHOOLEQ45009



Description :

My Block Mini Digital Dry Bath

Technical Specification :

My Block is the first and only mini heat block featuring a heated lid for protecting samples against unwanted condensation and evaporation.

Both the block and lid temperature can be selected independently in increments of 0.1°C, making this My Block the ideal choice for the most temperature sensitive applications, including immunoassays, lamp assays, enzyme reactions and other temperature hyper-sensitive procedures.

Ambient Range: 100°C	
Temp Accuracy:	

Temp. Increments:	
Temp. Uniformity:	
Dimensions: x 6 x 4.5 in./ 10.5 x 15 x 11.5cm (w x d x h)	
Block Dimensions: 4.7 x 7.1 x 3.2 cm	
Weight: 1 kg	
Electrical: Universal Voltage, 100 to 240V, 50-60Hz, 50W	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003