

Product Name :
Nauga Automatic Mounting Press

Product Code :
NLAB-ENGINEERINGLB84011



Description :

Nauga Automatic Mounting Press

Technical Specification :

Technical Specifications:

Pressure: 0~2MPa product pressure: 0~72MPa.

Shortest mounting time: 15 minutes.

Mould diameter (can be randomly selected 2 types): 25mm, 30mm, 40mm, 50mm.

Mounting temperature: 50-300°C.

Heating power: 1000W.

Voltage: 220V.

Weight: 120Kg.

Dimension: 630-520-520mm.

Features:

Microcomputer control, automatic Pre-heating, heating, applying pressure, cooling, etc.

Water-forced cooling, quickly get the sample.

Two demos can be done simultaneously.

Pressure, temperature, heating and cooling time can be set via touch screen.

Large size touch screen LCD, graphical HMI for friendly and easy operation.

Four kinds diameter moulds (25mm 30mm 40mm 50mm) can be chosen, easily change.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003