

Product Name :
Naugra Drill Tool Dynamometer

Product Code :
NLAB-ENGINEERINGLB15005



Description :

Naugra Drill Tool Dynamometer

Technical Specification :

It measures both the Thrust and Torque force of the drill bit and the torque produced on the work piece. It is used to establish drilling forces, study tool configuration and lubricant characteristics

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003