

Product Name :
Nauga Lab Abbe Refractometer

Product Code :
NLAB-LABORATORYEM75005



Description :

Nauga Lab Abbe Refractometer

Technical Specification :

Abbe Refractometer

Features:- of Abbe Refractometer

Measures dispersion

Horizontal loading of prism

Convenient reading angles

Used for measuring refractive index.

Range 1.3 to 1.7 with an accuracy of 0.001 direct on scale and 0.0001 by estimation.

Sugar percentage range of 0 to 95% with an accuracy of 1% on scale and 0.1 by estimation can be measured.

Complete with test piece, contact liquid and thermometer in wooden cabinet.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003